CUSTOM MOLDED EAR PLUGS - RADIANS CUSTOM MOLDED EAR PLUGS, BLUE

Do-It-Yourself Kit For A Custom Fit

Easy, do-it-yourself molded ear plugs provide a snug, permanent custom fit for hearing protection that blocks out harmful noise. Two-part, silicone kit mixes together to form a soft, pliable material used to mold to your ear canal. Dries in 10 minutes and cures in 3 hours. Made from safe, non-toxic and hypo-allergenic silicone, plus completely washable for long service life. Use alone or in combination with muffs for maximum protection.

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CUSTOM MOLDED EAR PLUGS SAFETY INSTRUCTIONS

Introduction

Thank you for choosing the CUSTOM MOLDED EAR PLUGS from RADIANS. These ear plugs are designed to provide effective hearing protection by blocking out harmful noise. This guide contains important safety instructions and information to ensure your safe and effective use of the product.

General Safety Guidelines

- Ensure you read and understand all instructions before use.
- Always follow the safety precautions outlined in this guide.
- Keep the product out of reach of children and vulnerable individuals.
- Check the product regularly for signs of wear or damage.
- Report any unsafe products or accidents to the relevant authorities.
- Stay informed about product recalls and safety updates via the EU's Safety Gate platform.

Specific Safety Precautions for Use

- **Age Restrictions:** These ear plugs are not recommended for children under 12 years of age without adult supervision.
- **Allergic Reactions:** While made from hypoallergenic silicone, discontinue use if you experience any irritation or allergic reactions.
- **Proper Fit:** Ensure that the ear plugs fit snugly in your ear canal to maximize noise reduction.
- **Combination Use:** For maximum protection, consider using the ear plugs in combination with earmuffs.
- **Moisture:** Avoid using the ear plugs in excessively wet environments to maintain their integrity.

Instructions for Installation and Usage

1. Preparation:

- Ensure your hands are clean and dry before handling the silicone kit.
- Gather all necessary materials from the kit.

2. Mixing the Silicone:

- Take the two parts of the silicone kit and mix them together thoroughly until a uniform color is achieved.

3. Molding the Ear Plugs:

- Immediately after mixing, take a small amount of the silicone and roll it into a ball.
- Insert the silicone ball into your ear canal and mold it to fit snugly. Ensure it is comfortable but secure.

4. Drying:

- Allow the ear plugs to dry for 10 minutes. Do not disturb them during this time.
- After 10 minutes, the ear plugs will be set but will continue to cure for a total of 3 hours. Avoid using them until fully cured.

5. Cleaning and Maintenance:

- The ear plugs are washable. Clean them with mild soap and water to maintain hygiene.
- Ensure they are completely dry before storing.

Disposal Instructions

- Dispose of the ear plugs responsibly. Do not throw them in regular household waste.
- Check local regulations for proper disposal methods for silicone products.

Contact Information for Further Support

For any questions or concerns regarding the CUSTOM MOLDED EAR PLUGS, please consult the manufacturer or your point of purchase for assistance.

Thank you for prioritizing your hearing safety with RADIANS CUSTOM MOLDED EAR PLUGS.

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