MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT BLACK

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050367
- Mfr. No.: HW-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.308kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650008

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Micro Trauma Kit NOW! Essentials Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This compact trauma kit is designed for everyday carry, ensuring you have essential lifesaving medical supplies readily available. This guide provides important safety instructions and general guidelines to ensure safe and effective use of the kit.

General Safety Guidelines

- Ensure that the Micro TKN is stored in a dry, cool place away from direct sunlight.
- Regularly check the contents of the kit to ensure that all supplies are intact and within their expiration dates.
- Keep the kit out of reach of children to prevent accidental ingestion or misuse of the medical supplies.
- In case of any allergic reactions to the supplies, discontinue use immediately and seek medical attention.
- Report any unsafe products or incidents related to the Micro TKN to the relevant authorities.

Specific Safety Precautions for Use

- Always read the instructions provided with each medical supply before use.
- Use trauma gloves when handling any blood or bodily fluids to prevent contamination.
- Ensure that the Quickclot dressing is applied according to the instructions for optimal bleeding control.
- When using the mini compression dressing, apply sufficient pressure to control bleeding without cutting off circulation.
- The hypothermia blanket should be used to maintain body temperature in cases of shock or hypothermia.
- Dispose of used medical supplies in accordance with local regulations regarding hazardous waste.

Instructions for Installation and Usage

Installation

1. Attaching the Kit:

- The Micro TKN can be mounted on a belt or MOLLE system.
- Ensure that the pouch is securely attached and does not move during activity.

Usage

1. Accessing Supplies:

- Use one hand or a single finger to pull the BLIP featured pull tabs to access supplies quickly.
- Practice accessing the kit in a safe environment to become familiar with its operation.

2. Using Medical Supplies:

- **Trauma Gloves:** Put on gloves before attending to any injuries.
- **Quickclot Dressing:** Open the package, apply directly to the wound, and hold in place to control bleeding.
- **Mini Compression Dressing:** Open the package and apply over the Quickclot dressing, securing it in place.
- **Hypothermia Blanket:** Unfold the blanket and wrap it around the individual in need, ensuring full coverage.
- **Flat Fold Tape:** Use to secure dressings in place.

Disposal Instructions

- Dispose of used medical supplies, including gloves, dressings, and tape, in accordance with local regulations.
- Do not throw away medical supplies in regular trash unless it is safe to do so.
- For expired or damaged supplies, follow local hazardous waste disposal guidelines.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW!, please refer to the product packaging for appropriate contact information. It is important to have access to support for any questions or concerns you may have regarding the use of this product.

By following these safety instructions, you can ensure the effective and safe use of the Micro Trauma Kit NOW! for yourself and others. Thank you for prioritizing safety and preparedness.

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