

MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MNT COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MNT COYOTE](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050368
- Mfr. No.: HW-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650015

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Safety Instructions](#)
- [About Us](#)

Micro Trauma Kit NOW! Safety Instructions

Introduction

Congratulations on your purchase of the Micro Trauma Kit NOW! (Micro TKN) by Blue Force Gear. This kit is designed for everyday carry and is equipped with essential lifesaving medical supplies. It is important to follow the safety instructions outlined in this guide to ensure safe and effective use of the kit.

General Safety Guidelines

- Ensure that all components of the kit are intact and in good condition before use.
- Keep the kit out of reach of children to prevent accidental use.
- Regularly check the expiration dates of medical supplies and replace any expired items immediately.
- Familiarize yourself with the contents of the kit and understand how to use each item before an emergency arises.
- Store the kit in a cool, dry place away from direct sunlight to preserve the integrity of the supplies.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Always wear gloves when providing first aid to protect yourself and the injured person.
 - Dispose of gloves properly after use to avoid contamination.
- Quickclot Bleeding Control Dressing:
 - Apply directly to the wound and maintain pressure to control bleeding.
 - Do not remove the dressing once applied, as this can disrupt clotting.
- Mini Compression Dressing:
 - Use as a secondary dressing for additional support on bleeding wounds.
 - Ensure the dressing is secure but not too tight, as this can restrict blood flow.
- Hypothermia Blanket:
 - Use the blanket to retain body heat in cases of shock or exposure.
 - Ensure the blanket is wrapped securely around the individual.
- Flat Fold Tape:
 - Use tape to secure dressings in place.
 - Ensure that the tape adheres well to avoid loosening during use.

Instructions for Installation and Usage

1. Deployment of the Kit:

- The Micro TKN can be deployed using one hand or a single finger by pulling the BLIP featured pull tabs.
- Ensure you are familiar with the deployment mechanism before an emergency situation arises.

2. Using the Medical Supplies:

- Open the kit and retrieve the necessary medical supplies based on the injury.
- Follow the specific instructions for each item as outlined above.
- If unsure about how to use any item, refer to additional resources or seek assistance from a trained professional.

3. After Use:

- Dispose of used medical supplies in accordance with local regulations.
- Restock any items that have been used or are nearing expiration to ensure the kit is always ready.

for use.

Disposal Instructions

- Dispose of expired or used medical supplies in accordance with local waste disposal regulations.
- Do not dispose of medical supplies in regular trash if they are biohazardous. Follow local guidelines for hazardous waste disposal.
- Ensure that the kit itself is stored properly or disposed of if damaged beyond repair.

Contact Information for Further Support

For any safety inquiries or concerns regarding the Micro Trauma Kit NOW!, please refer to the contact information provided by the retailer or manufacturer. It is important to have access to support for any questions or issues that may arise.

Conclusion

The Micro Trauma Kit NOW! is a vital tool for providing immediate first aid in emergencies. By following the safety instructions and guidelines outlined in this document, you can ensure the effective and safe use of the kit. Always stay prepared and informed to handle any medical situation that may arise.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk