MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MNT MULTICAM

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MNT MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050372
- Mfr. No.: BT-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650121

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Micro Trauma Kit NOW! Essentials Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Essentials Supplies. This kit is designed for everyday carry and is intended to provide essential lifesaving medical supplies. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the Micro Trauma Kit NOW! is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all supplies are intact and within their expiration dates.
- Keep the kit out of reach of children unless they are under adult supervision.
- In case of an emergency, always seek professional medical assistance after using the kit.
- Report any unsafe product experiences or accidents to the relevant authorities.
- Check for recall updates on the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to minimize the risk of infection.
- Use the Quickclot Bleeding Control Dressing only as directed to ensure proper application and effectiveness.
- Be cautious when using the Mini Compression Dressing; apply it firmly but do not overcompress, as this may restrict blood flow.
- The Hypothermia Blanket should be used to maintain body temperature; ensure it is wrapped securely around the individual in need.
- Flat Fold Tape should be used to secure dressings; ensure it adheres properly to avoid slippage.
- If you are unsure about any procedure, seek immediate medical assistance.

Instructions for Installation and Usage

1. Accessing the Kit:

- The Micro Trauma Kit NOW! can be accessed using the BLIP featured pull tabs. Pull from either the left or right side with one hand or a single finger.

2. Using the Supplies:

- **Trauma Gloves:**
 - Put on the gloves before handling any medical supplies or treating injuries.
- **Quickclot Bleeding Control Dressing:**
 - Apply directly to the wound and hold in place for at least three minutes to control bleeding.
- **Mini Compression Dressing:**
 - Place over the wound after applying Quickclot and secure with Flat Fold Tape.
- **Hypothermia Blanket:**
 - Wrap the individual securely, ensuring the reflective side is facing inward to retain body heat.
- **Flat Fold Tape:**
 - Use to secure dressings in place. Cut to the required length and apply firmly.

3. Repacking the Kit:

- After use, ensure all supplies are replaced and the kit is repacked for future use.

Disposal Instructions

- Dispose of used medical supplies, such as gloves and dressings, in accordance with local regulations for medical waste.
- The kit itself should be disposed of responsibly if damaged or beyond repair. Check with local waste management services for proper disposal methods.

Contact Information for Further Support

For any inquiries or concerns regarding the Micro Trauma Kit NOW! Essentials Supplies, please refer to the EU contact point for safety inquiries.

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