MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT COYOTE

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050374
- Mfr. No.: HW-TKN-MTKN-PRO
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650060

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Micro Trauma Kit NOW! Pro Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies. This guide provides essential safety instructions and guidelines to ensure the safe and effective use of your trauma kit. Please read this manual carefully before using the product.

General Safety Guidelines

- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all supplies are intact and within their expiration dates.
- Keep the kit out of reach of children unless under adult supervision.
- Use the kit only for its intended purpose, which is to provide basic first aid in emergency situations.
- Familiarize yourself with the contents and their uses before an emergency arises.
- Report any unsafe products or accidents to the appropriate authorities.
- Stay updated on recall notices and safety alerts through the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to prevent contamination.
- Use the Quickclot Combat Gauze LE only as directed in emergency situations.
- Ensure that the SWATT Pressure Device is applied correctly to avoid further injury.
- The Hypothermia Blanket should be used to maintain body temperature; do not use if the patient is already overheating.
- Be cautious when using the kit in extreme weather conditions that may affect the integrity of the supplies.

Instructions for Installation and Usage

1. **Familiarization:** Before use, open the pouch and familiarize yourself with each item included in the kit.
2. **Deployment:**
 - The Micro TKN can be deployed using one hand. Pull the BLIP featured pull tabs from either the left or right side.
3. **Using Medical Supplies:**
 - **Trauma Gloves:** Put on the gloves before attending to any injuries.
 - **Quickclot Combat Gauze LE:** Apply directly to the wound and maintain pressure.
 - **SWATT Pressure Device:** Wrap around the injury site and secure it tightly.
 - **Hypothermia Blanket:** Wrap the patient to prevent heat loss.
4. **PostUse:** After using any supplies, dispose of them properly and restock the kit as needed.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Do not throw away the kit itself; instead, ensure it is restocked and ready for future use.
- Check for any specific disposal instructions on individual items within the kit.

Contact Information for Further Support

For any inquiries regarding safety or product concerns, please reach out to the designated EU contact point associated with this product.

Thank you for prioritizing safety and preparedness with the Micro Trauma Kit NOW! Pro Supplies. Your attention to these guidelines will help ensure effective use in critical situations.

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