

# MICRO TRAUMA KIT NOW! PRO SUPPLIES - BLUE FORCE GEAR

## MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components - the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)



### Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT COYOTE](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050374
- Mfr. No.:
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650060

### Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Pro Supplies Safety Instruction Guide](#)
- [About Us](#)

# Micro Trauma Kit NOW! Pro Supplies Safety Instruction Guide

## Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies. This guide provides essential safety instructions to ensure safe use and handling of the kit. The Micro Trauma Kit NOW! is designed for everyday carry and is suitable for law enforcement professionals, prepared citizens, and hunters. Please read this guide carefully to understand the safety measures and proper usage of the kit.

## General Safety Guidelines

To ensure your safety and the safety of others, please follow these guidelines:

- Always inspect the kit for any visible damage or missing items before use.
- Keep the kit out of reach of children and vulnerable individuals when not in use.
- Use the kit only for its intended purpose, which is to provide first aid in emergencies.
- Ensure that all components are stored in their designated places within the kit.
- Regularly check and replace expired or used items in the kit.
- Report any unsafe products or incidents to the appropriate authorities.
- Stay informed about product recalls by checking the EU Safety Gate platform.

## Specific Safety Precautions for Use

When using the Micro Trauma Kit NOW!, please observe the following precautions:

- **Trauma Gloves:** Always wear trauma gloves when providing first aid to prevent contamination and protect yourself from bodily fluids.
- **Quickclot Combat Gauze:** Use only as directed. Do not apply directly to open wounds without first following proper first aid procedures.
- **SWATT Pressure Device:** Ensure the device is applied correctly to control bleeding. If unsure, seek assistance from a trained professional.
- **Hypothermia Blanket:** Use the blanket to maintain body temperature in case of exposure to cold; ensure it is wrapped securely around the individual.
- **OneHanded Operation:** Practice using the BLIP pull tabs with one hand to ensure proficiency in emergency situations.

## Instructions for Installation and Usage

Follow these instructions for the proper installation and usage of the Micro Trauma Kit NOW!:

### 1. Mounting the Kit:

- Attach the MOLLE or beltmounted pouch securely to your belt or gear using the TenSpeed technology.
- Ensure the pouch is firmly in place and does not shift during movement.

### 2. Accessing Supplies:

- To access the supplies, pull the BLIP featured pull tabs located on either side of the pouch.
- Practice pulling the tabs with one hand to ensure you can access supplies quickly in an emergency.

### 3. Using the Supplies:

- Follow the specific instructions for each component as outlined in the safety precautions above.
- In case of serious injuries, do not hesitate to call emergency services for professional assistance.

#### **4. After Use:**

- Dispose of used items, such as gloves and gauze, in accordance with local regulations.
- Restock any used or expired items promptly to ensure the kit is always ready for use.

## **Disposal Instructions**

When disposing of the Micro Trauma Kit NOW! or its components, please adhere to the following guidelines:

- Dispose of used medical supplies in accordance with local health regulations.
- Recycle packaging materials where possible.
- Do not dispose of the kit in regular household waste if it contains hazardous materials. Check local regulations for proper disposal methods.

## **Contact Information for Further Support**

For any questions or concerns regarding the Micro Trauma Kit NOW! Pro Supplies, please refer to the manufacturer's contact information provided with the product packaging. It is important to have access to support for any safety inquiries or product issues.

Thank you for prioritizing safety with the Micro Trauma Kit NOW! Pro Supplies. Your preparedness can make a significant difference in emergency situations.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)