MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT ODG

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT ODG
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100743
- Mfr. No.: HW-TKN-MTKN-ESS-OD
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073650039

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Micro Trauma Kit NOW! Essentials Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This kit is designed to provide essential lifesaving medical supplies for everyday carry. Whether you are a law enforcement professional, a prepared citizen, or a hunter, this guide will help you understand how to use the kit safely and effectively.

General Safety Guidelines

- Ensure that all components of the Micro TKN are intact and in good condition before use.
- Keep the kit out of reach of children and vulnerable individuals.
- Regularly check for updates on product recalls or safety alerts via the EU Safety Gate platform.
- Report any unsafe products or incidents to local authorities.
- Always use the kit in accordance with the instructions provided.

Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear gloves when handling blood or bodily fluids to prevent infection.
- **Bleeding Control Dressing:** Apply pressure directly to the wound and secure the dressing firmly. If bleeding does not stop, seek professional medical help immediately.
- **Mini Compression Dressing:** Use only as instructed. If the dressing becomes soaked with blood, replace it promptly.
- **Hypothermia Blanket:** Use the blanket to retain body heat in cold environments. Ensure it is wrapped securely around the individual.
- **Flat Fold Tape:** Use tape to secure dressings in place. Avoid applying tape directly to the skin to prevent irritation.

Instructions for Installation and Usage

1. Accessing the Kit:

- The Micro TKN is designed for quick access. Use the BLIP pull tabs to open the pouch with one hand or a single finger.

2. Using the Components:

- **Trauma Gloves:** Remove from the pouch and put on before treating any injuries.
- **Quickclot Bleeding Control Dressing:**
 - Open the packaging and apply directly to the wound.
 - Hold in place for at least three minutes to allow clotting.
- **Mini Compression Dressing:**
 - Place over the Quickclot dressing if additional compression is needed.
 - Secure with Flat Fold Tape.
- **Hypothermia Blanket:**
 - Unfold the blanket and wrap it around the individual.
 - Ensure the shiny side is facing the individual for maximum heat retention.
- **Flat Fold Tape:**
 - Cut the desired length and apply to secure dressings in place.

3. PostUse:

- After use, dispose of all singleuse items, including gloves and dressings, in accordance with local regulations.
- Replenish the kit with necessary supplies as soon as possible.

Disposal Instructions

- Dispose of used medical supplies in accordance with local health regulations.
- Do not dispose of any part of the kit in regular household waste if it contains biohazardous materials.
- Recycle any packaging materials where possible.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW!, please refer to the contact information provided in the product packaging or user manual. Ensure that you have the product details ready when seeking assistance.

Thank you for prioritizing safety with the Micro Trauma Kit NOW!

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