MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MOUNT ODG

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES BELT MOUNT ODG
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100744
- Mfr. No.: BT-TKN-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073650138

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Micro Trauma Kit NOW! Essentials Supplies Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This compact trauma kit is designed for everyday carry, providing essential lifesaving medical supplies on your belt. The Micro TKN is ideal for law enforcement professionals, prepared citizens, and hunters. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the Micro TKN is stored in a location that is easily accessible.
- Always check the contents of the kit regularly to ensure all supplies are intact and within their expiration dates.
- Use the kit only for its intended purpose: to provide first aid in emergency situations.
- Keep the kit out of reach of children unless supervised by an adult.
- Report any unsafe products or incidents to the relevant authorities.

Specific Safety Precautions for Use

- Always wear gloves when handling blood or bodily fluids to minimize the risk of infection.
- Be aware of your surroundings when providing first aid; ensure your safety and that of the injured person.
- Do not use any supplies that appear damaged or expired.
- Follow the instructions provided with each medical supply for proper use.
- If you are unsure about how to use any item in the kit, seek guidance from a qualified medical professional.

Instructions for Installation and Usage

1. Mounting the Kit:

- Attach the Micro TKN pouch to your belt using the MOLLE or Belt mount system.
- Ensure that the pouch is securely fastened and does not obstruct movement.

2. Accessing Supplies:

- To access supplies, pull the BLIP featured pull tabs with one hand or a single finger from either side of the pouch.
- Practice accessing the kit in a safe environment to become familiar with the process.

3. Using the Medical Supplies:

- **Trauma Gloves:** Always wear gloves before providing first aid. Dispose of gloves properly after use.
- **Quickclot Bleeding Control Dressing:** Apply directly to the wound to control bleeding. Follow the instructions on the packaging for best results.
- **Mini Compression Dressing:** Use to provide compression on wounds. Ensure it is applied firmly but not so tight as to cut off circulation.
- **Hypothermia Blanket:** Use to prevent heat loss in a person who is injured or in shock. Wrap the individual snugly in the blanket.
- **Flat Fold Tape:** Use to secure dressings in place. Ensure the area is clean and dry before application.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Do not throw away used gloves or dressings in regular trash; use designated biohazard disposal methods if available.
- Check the expiration dates of supplies regularly and dispose of any expired items appropriately.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW!, please refer to the manufacturer or retailer where you purchased the product for assistance. Always ensure that you have access to a reliable source for safety inquiries.

Thank you for your attention to these safety instructions. Your safety and the safety of others depend on the proper use of this kit.

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