MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - BLACK

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal.

Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101011
- Mfr. No.: HW-TKN-PLUS-ADV-BK
- Size: 1 Person
- Delivery weight: 0.227kg
- Shipping height: 114mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653115

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Micro Trauma Kit NOW! PLUS+ Safety Instruction Guide

Introduction

Thank you for purchasing the Micro Trauma Kit NOW! PLUS+. This advanced medical pouch is designed to provide essential first aid supplies in a compact and accessible format. This guide is intended to ensure your safety and the effective use of the product in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- **Product Safety:** Ensure safe use of the Micro Trauma Kit by following the instructions provided. This kit is designed for nonfood use and is compliant with EU safety standards.
- **Enhanced Recalls:** Stay informed about any recalls. In the event of a recall, standardized notice will be provided, and free remedies will be available for dangerous products.
- **Online Shopping:** If purchased online, ensure the platform meets safety requirements similar to physical stores.
- **Special Consumer Focus:** Extra care should be taken when using this kit around vulnerable groups, such as children.
- **EU Contact Point:** For any safety inquiries, refer to the designated EU contact point for this product.
- **Rapid Alerts:** Be aware of updates regarding unsafe products through the EU's Safety Gate system.

Specific Safety Precautions for Use

- **Read Instructions:** Always read and understand the instructions before using the Micro Trauma Kit.
- **Check Contents:** Regularly check the contents of the kit to ensure all items are present and not expired.
- **Use Gloves:** Always wear trauma gloves when administering first aid to prevent contamination and protect yourself.
- **Proper Disposal:** Dispose of used medical supplies according to local regulations to prevent health hazards.
- **Age Restrictions:** This kit is not intended for use by children without adult supervision.
- **Storage Conditions:** Store the kit in a cool, dry place away from direct sunlight and moisture.
- **Emergency Situations:** Use the kit only in emergency situations and follow the instructions for each item carefully.

Instructions for Installation and Usage

Installation

1. **MOLLE Attachment:**
 - Ensure the MOLLE med kit mount is compatible with your equipment (e.g., plate carriers, vests, battle belts).
 - Securely attach the kit using the MOLLE system by threading the straps through the webbing.

Usage

1. Accessing Supplies:

- Open the pouch by unzipping or unclipping the closures.
- Identify the required supplies based on the nature of the injury or emergency.

2. Using Trauma Gloves:

- Roll out the trauma gloves and wear them before handling any medical supplies.
- Dispose of gloves properly after use.

3. Applying QuikClot Dressing:

- Open the QuikClot Bleeding Control Dressing and apply it directly to the wound.
- Follow the instructions on the packaging for proper application.

4. Using Emergency Trauma Dressing:

- Remove the trauma dressing from the pouch.
- Apply it to the affected area, ensuring adequate coverage.

5. Using Hypothermia Blanket:

- Unwrap the hypothermia blanket and place it over the patient to maintain body temperature.

6. Using Flat Fold Tape:

- Cut the desired length of tape and use it to secure dressings in place.

7. Using Pressure Device:

- Apply the pressure device to control bleeding as instructed.

8. Using Chest Seal:

- Open the chest seal package and apply it to the chest wound as directed.

9. Using Nasopharyngeal Airway:

- Insert the nasopharyngeal airway as instructed to maintain an open airway.

10. Using Decompression Needle:

- Use the decompression needle only if you are trained to do so, following the specific instructions.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local regulations.
- Ensure that sharps (e.g., needles) are disposed of in a designated sharps container.
- Do not throw away medical supplies in regular trash to avoid potential health risks.

Contact Information for Further Support

For any safety inquiries or additional support regarding the Micro Trauma Kit NOW! PLUS+, please refer to the designated EU contact point for this product.

By following these guidelines, you can ensure the safe and effective use of the Micro Trauma Kit NOW! PLUS+. Thank you for prioritizing safety and preparedness.

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