MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - COYOTE

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal.

Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101012
- Mfr. No.: HW-TKN-PLUS-ADV-CB
- Size: 1 Person
- Delivery weight: 0.227kg
- Shipping height: 305mm
- Shipping width: 83mm
- Shipping length: 254mm
- UPC: 810073653122

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Micro Trauma Kit NOW! PLUS+ Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! PLUS+. This advanced medical pouch is designed for safety and efficiency in emergency situations. Please read this safety instruction guide carefully to ensure safe use and compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure that the Micro Trauma Kit is kept out of reach of children and vulnerable individuals to prevent accidental use.
- Regularly check the kit for any expired or damaged items. Replace as necessary to maintain effectiveness.
- Use the kit only for its intended purpose. It is designed for emergency medical situations and should not be used as a substitute for professional medical assistance.
- Familiarize yourself with the contents of the kit and how to use each item effectively.
- Report any unsafe products or incidents to the relevant authorities and check for recall updates via the EU's Safety Gate platform.

Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear gloves when providing first aid to prevent contamination and protect both the caregiver and the patient.
- **QuikClot Products:** Follow the instructions provided with the QuikClot products carefully to ensure proper bleeding control.
- **Emergency Trauma Dressing:** Ensure that the dressing is applied correctly to manage wounds effectively. If unsure, seek professional assistance.
- **Hypothermia Blanket:** Use the blanket to prevent heat loss in patients who are experiencing shock or hypothermia.
- **Nasopharyngeal Airway:** Only use if trained in its application. Incorrect use can cause harm.
- **Decompression Needle:** This should only be used by trained individuals. Improper use can lead to severe injury.

Instructions for Installation and Usage

1. **Familiarization:** Open the Micro Trauma Kit and review the contents. Ensure you understand the purpose and usage of each item.
2. **Carrying the Kit:** Attach the kit to a MOLLE compatible surface, such as a plate carrier or vest, ensuring it is secure for easy access.
3. **Using the Kit:**
 - In an emergency, quickly access the kit and remove the necessary items.
 - Follow the specific instructions for each item as provided in the kit or on the packaging.
 - Always prioritize calling for professional medical assistance in serious situations.
4. **PostUse:** After using any items, replace them immediately to ensure the kit is ready for future use. Dispose of any used materials according to local regulations.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations. Many items may be considered biohazardous waste.
- Check with local waste management authorities for proper disposal methods for medical waste.
- Ensure that the Micro Trauma Kit itself is kept intact and does not become damaged during disposal.

Contact Information for Further Support

For safety inquiries or support regarding the Micro Trauma Kit NOW! PLUS+, please refer to the contact information provided on the product packaging or visit the manufacturer's website.

Thank you for prioritizing safety and preparedness with the Micro Trauma Kit NOW! PLUS+. Stay safe and be prepared for any emergency situation.

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