MICRO TRAUMA KIT NOW! - PLUS+ - MICRO TRAUMA KIT NOW!-PLUS+ - MOLLE-ADVANCED-RANGER GREEN

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-PLUS+ - MOLLE-ADVANCED-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101014
- Mfr. No.: HW-TKN-PLUS-ADV-RG
- Size: 1 Person
- Delivery weight: 0.526kg
- Shipping height: 102mm
- Shipping width: 241mm
- Shipping length: 279mm
- UPC: 810073653146

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MICRO TRAUMA KIT NOW! PLUS+ Safety Instructions

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+. This comprehensive medical pouch is designed to provide essential first aid supplies for emergencies. To ensure your safety and the effective use of this product, please read the following safety instructions carefully.

General Safety Guidelines

- Ensure the product is used only for its intended purpose as a first aid kit.
- Keep the kit out of reach of children to prevent misuse or accidental ingestion of components.
- Regularly check the contents of the kit to ensure all items are present and within their expiration dates.
- Store the kit in a cool, dry place away from direct sunlight and extreme temperatures.
- Familiarize yourself with the contents of the kit and their uses before an emergency arises.
- Report any unsafe products or accidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear Trauma Gloves when handling blood or bodily fluids to reduce the risk of infection.
- Use the QuikClot Bleeding Control Dressing only as instructed on the packaging.
- Ensure that the Emergency Trauma Dressing is applied correctly to control bleeding effectively.
- The Hypothermia Blanket should be used to maintain body temperature in cases of shock or exposure.
- Be cautious when using the Decompression Needle; improper use can cause injury.
- The Nasopharyngeal Airway should only be used by trained individuals to ensure proper insertion.

Instructions for Installation and Usage

Installation

- The MICRO TRAUMA KIT NOW! PLUS+ can be mounted on MOLLEcompatible surfaces, such as plate carriers, vests, and battle belts.
- Ensure that the MOLLE straps are securely fastened to prevent the kit from detaching during use.

Usage

1. **Assess the Situation:** Before using the kit, evaluate the scene for safety and the nature of injuries.
2. **Wear Gloves:** Always put on Trauma Gloves before providing first aid.
3. **Control Bleeding:** Use the QuikClot Bleeding Control Dressing or Emergency Trauma Dressing as needed.
4. **Maintain Body Temperature:** If the individual is experiencing shock, wrap them in the Hypothermia Blanket.
5. **Airway Management:** If necessary, use the Nasopharyngeal Airway to maintain an open airway.
6. **Seek Professional Help:** After administering first aid, ensure that professional medical assistance is sought.

Disposal Instructions

- Dispose of used medical supplies, such as gloves and dressings, in accordance with local regulations.
- Do not dispose of the kit in regular household waste if it contains hazardous materials.
- Check local guidelines for the proper disposal of medical waste.

Contact Information for Further Support

For safety inquiries or to report unsafe products, please refer to the official EU Safety Gate platform or your local consumer safety authority.

By following these safety instructions, you can ensure the effective use of the MICRO TRAUMA KIT NOW! PLUS+ and contribute to a safer environment for yourself and others.

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