

MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW!-PLUS+ - MOLLE-ADVANCED-RANGER GREEN

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW!-PLUS+ - MOLLE-ADVANCED-RANGER GREEN](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101014
- Mfr. No.: HW-TKN-PLUS-ADV-RG
- Size: 1 Person
- Delivery weight: 0.526kg
- Shipping height: 102mm
- Shipping width: 241mm
- Shipping length: 279mm
- UPC: 810073653146

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Safety Instruction Guide for MICRO TRAUMA KIT NOW! PLUS+](#)
- [About Us](#)

Safety Instruction Guide for MICRO TRAUMA KIT NOW! PLUS+

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+ by BLUE FORCE GEAR. This guide provides essential safety instructions to ensure the safe and effective use of your trauma kit. It is important to read and understand these instructions to maximize your safety and the safety of others.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place, away from direct sunlight and moisture.
- Regularly check the contents of the kit for expiration dates and signs of damage.
- Keep the kit out of reach of children and vulnerable individuals unless supervised.
- Familiarize yourself with the contents and their uses prior to an emergency situation.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling medical supplies to prevent contamination.
- Use QuikClot products according to the instructions provided in the packaging.
- Ensure all dressings are sterile before application.
- Do not use any items in the kit that appear damaged or expired.
- Be cautious when using sharp objects, such as scissors or needles, to avoid injury.
- If you are unsure about how to use any item in the kit, seek professional assistance.

Instructions for Installation and Usage

1. Preparation

- Open the trauma kit and lay out all items on a clean surface.
- Identify the items you will need based on the situation.

2. Using the Trauma Gloves

- Roll out one pair of trauma gloves from the kit.
- Put on the gloves before handling any medical supplies.

3. Applying Dressings

- For bleeding control, use the QuikClot Bleeding Control Dressing:
 - Open the package carefully.
 - Apply the dressing directly to the wound and apply pressure.
- For larger wounds, use the Emergency Trauma Dressing:
 - Open the dressing and place it over the wound.
 - Secure it in place using the attached ties.

4. Using the Chest Seal

- In case of a chest injury, use the Full Size Chest Seal:
 - Ensure the area around the wound is clean.
 - Peel off the backing and apply the seal over the wound, ensuring a tight fit.

5. Managing Hypothermia

- If hypothermia is suspected, use the Hypothermia Blanket:
 - Unroll the blanket and wrap it around the individual to provide warmth.

6. Using the Nasopharyngeal Airway

- If airway obstruction is suspected, use the Nasopharyngeal Airway:
 - Lubricate the airway if possible.
 - Gently insert the airway into the nostril until resistance is met.

7. Using the Decompression Needle

- In cases of tension pneumothorax, use the Decompression Needle:
 - Identify the correct site on the chest.
 - Insert the needle quickly and with caution to relieve pressure.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not dispose of any items in regular household waste if they are contaminated or hazardous.
- Always check local guidelines for the disposal of sharps and medical waste.

Contact Information for Further Support

For any questions or concerns regarding the MICRO TRAUMA KIT NOW! PLUS+, please refer to the product packaging for the appropriate contact information.

By following these safety instructions, you can ensure the effective use of your MICRO TRAUMA KIT NOW! PLUS+ and contribute to a safer environment for yourself and those around you. Thank you for your attention to these important guidelines.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk