TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - BLACK

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101016
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.186kg
- Shipping height: 114mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653214

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Deutsch: TRAUMA KIT NOW! SMALL BLUE FORCE GEAR TRAUMA KIT NOW! SMALL MOLLE ESSENTIALS SUPPLIES BLACK Sicherheitshinweise](#)
- [English: Trauma Kit NOW! SMALL Safety Instruction Guide](#)
- [Français: Guide de Sécurité pour la Trousse de Secours TRAUMA KIT NOW! SMALL](#)
- [Italiano: Guida alle Istruzioni di Sicurezza per il Trauma Kit NOW! SMALL](#)
- [Svenska: Säkerhetsinstruktioner för TRAUMA KIT NOW! SMALL](#)
- [Český: Návod k použití: TRAUMA KIT NOW! SMALL](#)

TRAUMA KIT NOW! SMALL BLUE FORCE GEAR

TRAUMA KIT NOW! SMALL MOLLE ESSENTIALS

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Trauma Kit NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL by Blue Force Gear. This compact medical pouch is designed to provide essential supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that all components of the Trauma Kit NOW! SMALL are intact and undamaged before use.
- Keep the kit out of reach of children unless under adult supervision.
- Store the kit in a cool, dry place to maintain the integrity of the medical supplies.
- Regularly check the expiration dates of the medical supplies and replace any expired items.
- Familiarize yourself with the contents of the kit and how to use each item before an emergency occurs.
- In case of a medical emergency, always seek professional medical assistance as soon as possible.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Use gloves to protect yourself and the injured person from infection.
 - Dispose of gloves properly after use to prevent contamination.
- QuikClot Bleeding Control Dressing:
 - Apply directly to the wound and apply pressure as needed.
 - Do not remove the dressing once applied; seek medical assistance.
- Emergency Trauma Dressing:
 - Use as a primary dressing for severe wounds.
 - Follow the instructions on the packaging for proper application.
- Hypothermia Blanket:
 - Use to maintain body temperature in case of shock or exposure.
- Flat Fold Tape:
 - Use to secure dressings in place.
 - Ensure the skin is clean and dry before application.
- Compressed Gauze:
 - Use for wound packing to control bleeding.
 - Apply pressure to the wound site after packing.
- Full Size Chest Seal:
 - Use for open chest wounds to prevent air from entering the chest cavity.
 - Follow the instructions on the packaging for proper application.
- Nasopharyngeal Airway:
 - Only trained individuals should use this device.
 - Ensure proper sizing and lubrication before use.
- Decompression Needle:
 - Only trained individuals should use this device.
 - Follow proper medical protocols for use.

Instructions for Installation and Usage

1. **Locate the Trauma Kit NOW! SMALL** in your designated storage area.
2. **Open the kit** and familiarize yourself with the contents.
3. **In an emergency**, assess the situation and determine the necessary supplies needed from the kit.
4. **Follow the specific safety precautions** listed above for each item you intend to use.
5. **After use**, ensure all items are properly disposed of or cleaned, as applicable.
6. **Restock the kit** with any used or expired items to ensure readiness for future emergencies.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not throw away used gloves, dressings, or other medical waste in regular trash.
- Check with local health authorities for proper disposal methods of hazardous medical waste.

Contact Information for Further Support

For any inquiries regarding the Trauma Kit NOW! SMALL or to report unsafe products, please refer to the EU Safety Gate platform for updates and further assistance.

Thank you for prioritizing safety with the Trauma Kit NOW! SMALL. Your preparedness can make a significant difference in emergency situations.

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