TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - COYOTE

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101017
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 229mm
- Shipping length: 254mm
- UPC: 810073653221

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TRAUMA KIT NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This comprehensive medical kit is designed to provide essential supplies for trauma care in a compact and portable format. Please read this safety instruction guide carefully to ensure safe and effective use of the kit.

General Safety Guidelines

- **Product Safety:** This product is designed for safe use and complies with EU safety regulations. Always follow the instructions provided to minimize risks.
- **Enhanced Recalls:** In the event of a product recall, you will be notified through standardized recall notices. Please follow the instructions provided for any remedies.
- **Online Shopping Safety:** If purchased online, ensure that the platform complies with safety requirements similar to physical stores.
- **Special Consumer Focus:** This kit is not intended for children. Keep out of reach of children to prevent accidental use or injury.
- **EU Contact Point:** For any safety inquiries, please refer to the contact information provided with your product.
- **Rapid Alerts:** Stay informed about product safety updates through the EU's Safety Gate platform. Report any unsafe products or accidents to local authorities.

Specific Safety Precautions for Use

- **Check Contents:** Before use, verify that all components of the kit are present and in good condition.
- **Use as Intended:** This kit is designed for trauma care. Do not use for any other purpose.
- **Follow Instructions:** Always follow the instructions included for each item within the kit to ensure proper use.
- **Hygiene Practices:** Use gloves provided in the kit when treating wounds to minimize infection risk.
- **Avoid Direct Sunlight:** Store the kit in a cool, dry place away from direct sunlight to maintain the integrity of the supplies.
- **Monitor for Expiration:** Regularly check the expiration dates of items in the kit and replace as necessary.

Instructions for Installation and Usage

- **Opening the Kit:**
 - Use the MOLLE attachment system to securely attach the kit to your gear or bag.
 - Ensure that the kit is easily accessible in case of emergencies.
- **Using the Supplies:**
 - **Trauma Gloves:** Wear gloves when handling open wounds or bodily fluids to prevent contamination.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound and apply pressure to control bleeding.
 - **Emergency Trauma Dressing:** Place over the wound and secure in place using the attached ties.
 - **Hypothermia Blanket:** Wrap the patient to maintain body temperature in case of shock.
 - **Flat Fold Tape:** Use to secure dressings in place or to close small wounds.
 - **Compressed Gauze:** Apply to wounds to absorb blood and promote clotting.
 - **Full Size Chest Seal:** For chest injuries, apply the seal over the wound to prevent air from entering the chest cavity.
 - **Nasopharyngeal Airway:** Insert as needed to maintain an open airway in unresponsive patients.
 - **Decompression Needle:** Use only if trained, for severe chest injuries with signs of tension pneumothorax.

Disposal Instructions

- **Expired or Used Items:** Dispose of any expired or used medical supplies in accordance with local regulations for medical waste.
- **Packaging:** Recycle the outer packaging where possible, following local recycling guidelines.

Contact Information for Further Support

For any questions or concerns regarding the use of the TRAUMA KIT NOW! SMALL, please refer to the contact information provided with your product packaging.

Thank you for your attention to these safety instructions. Your safety and wellbeing are our top priority.

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