TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW!

SMALL-MOLLE-ESSENTIALS SUPPLIES-RANGER GREEN

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101019
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 51mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653245

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Safety Instruction Guide for TRAUMA KIT NOW! SMALL

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This medical pouch is designed to provide essential first aid supplies in a compact format. It is important to read and understand these safety instructions to ensure safe and effective use of the kit.

General Safety Guidelines

- **Product Safety:** This product is designed for safe use and complies with EU safety standards. Ensure that you follow all safety guidelines to minimize risks.
- **Enhanced Recalls:** In the event of a safety recall, you will be notified through standardized recall notices. Please follow the instructions provided for returning or exchanging the product.
- **Online Shopping:** If you purchased this kit online, you are entitled to the same safety protections as if you bought it in a physical store.
- **Special Consumer Focus:** Extra care should be taken when using this kit around vulnerable groups, such as children. Always supervise their use of medical supplies.
- **EU Contact Point:** For any safety inquiries, please consult the relevant EU-based contact point associated with this product.
- **Rapid Alerts:** Stay informed about product safety updates through the EU's Safety Gate platform.

Specific Safety Precautions for Use

- **Keep out of reach of children:** This kit contains items that may pose a choking hazard or injury risk.
- **Check contents regularly:** Ensure that all items are present and in good condition before use. Replace any expired or damaged items.
- **Use only as directed:** Follow the instructions for each item included in the kit. Misuse can lead to ineffective treatment or injury.
- **Store in a cool, dry place:** Keep the kit away from excessive heat or moisture to maintain the integrity of the supplies.
- **Be aware of allergies:** Check for any personal allergies to materials included in the kit, such as latex in gloves.

Instructions for Installation and Usage

- **Opening the Kit:**
 - Use the single-motion pullout design to access the contents quickly.
- **Utilizing the Supplies:**
 - **Trauma Gloves:** Always wear gloves when treating wounds to minimize contamination.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound to control bleeding. Follow the instructions on the packaging.
 - **Emergency Trauma Dressing:** Use in the event of severe bleeding. Place over the wound and secure in place.
 - **Hypothermia Blanket:** Wrap the injured person to retain body heat in case of shock or exposure.
 - **Flat Fold Tape:** Use to secure dressings or gauze in place.
 - **Compressed Gauze:** Use to pack wounds as necessary.
 - **Full Size Chest Seal:** Apply to open chest wounds to prevent air from entering the chest cavity.
 - **Nasopharyngeal Airway:** Use to maintain an open airway in unconscious patients.
 - **Decompression Needle:** Use only if trained to do so. This is a critical procedure for tension pneumothorax.

Disposal Instructions

- **Dispose of used items properly:** Follow local regulations for disposing of medical waste.
- **Check expiration dates:** Replace expired items in the kit and dispose of them according to local guidelines.
- **Do not reuse singleuse items:** Items such as gloves, dressings, and gauze should only be used once.

Contact Information for Further Support

For any questions or concerns regarding the TRAUMA KIT NOW! SMALL, please refer to the manufacturer's website or the retailer from whom you purchased the product for further assistance.

By following these safety instructions, you can ensure the effective and safe use of your TRAUMA KIT NOW! SMALL. Thank you for prioritizing safety and preparedness.

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