TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW!

SMALL-MOLLE-ESSENTIALS SUPPLIES-WOLF GRAY

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101020
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.442kg
- Shipping height: 102mm
- Shipping width: 178mm
- Shipping length: 254mm
- UPC: 810073653252

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Trauma Kit NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL from Blue Force Gear. This compact medical kit is designed for quick access to essential trauma supplies. It is crucial to understand how to use this kit safely and effectively. This guide provides important safety instructions and guidelines to ensure the safe use of your trauma kit.

General Safety Guidelines

- Always inspect the kit before use to ensure all components are intact and undamaged.
- Keep the kit out of reach of children unless supervised by an adult.
- Store the kit in a cool, dry place to prevent damage to the contents.
- Do not use any items in the kit after their expiration date.
- Familiarize yourself with the contents and their intended use before an emergency occurs.
- In case of a medical emergency, seek professional medical assistance immediately.

Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear gloves when treating wounds to prevent infection. Dispose of gloves properly after use.
- **QuikClot Products:** Follow instructions carefully when using QuikClot products for bleeding control. Ensure the area is clean before application.
- **Emergency Trauma Dressing:** Use the dressing only as directed. Do not apply excessive pressure; allow for blood flow if needed.
- **Hypothermia Blanket:** Use the blanket to retain body heat in cases of shock or hypothermia. Ensure it is wrapped properly around the person in need.
- **Compressed Gauze:** Apply gauze directly to the wound and hold in place. Change as necessary to control bleeding.
- **Chest Seal:** Ensure the chest seal is applied correctly to prevent air from entering the chest cavity in case of a penetrating chest wound.
- **Nasopharyngeal Airway:** Use only if trained in its application. Improper use can cause injury.
- **Decompression Needle:** This should only be used by trained personnel. Incorrect use can lead to serious injury.

Instructions for Installation and Usage

1. **Familiarization:** Before an emergency arises, take time to familiarize yourself with the kit's contents and their uses.
2. **Accessing the Kit:** Open the MOLLE pouch using the single-motion pullout design for quick access to supplies.
3. **Using Supplies:**
 - For each item, follow the specific usage instructions provided with the product.
 - In case of severe bleeding, prioritize the use of QuikClot and the Emergency Trauma Dressing.
4. **After Use:** Clean and disinfect any reusable supplies after use. Replace any used items promptly to maintain readiness.
5. **Regular Checks:** Periodically check the kit to ensure all items are present and within their expiration dates. Replace as necessary.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations.
- For medical waste, use designated disposal methods to ensure safety and compliance with health regulations.
- Recycle packaging materials where possible.

Contact Information for Further Support

For any inquiries regarding safety or product use, please refer to the manufacturer's contact details provided with your product. Ensure you have the product details ready for quicker assistance.

Conclusion

Your safety and the effective use of the Trauma Kit NOW! SMALL are our top priorities. Following these guidelines will help ensure that you are prepared to handle emergencies safely and effectively. Always prioritize professional medical help in any emergency situation.

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