TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES - BLACK

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES - BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101021
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.227kg
- Shipping height: 102mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653269

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TRAUMA KIT NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This compact trauma kit is designed to provide essential medical supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace expired items.
- Keep the kit out of reach of children unless supervised by an adult.
- If you experience any allergic reactions to items in the kit, seek medical attention immediately.
- Report any unsafe products or accidents to the appropriate authorities.
- Stay informed about product recalls through the EU Safety Gate platform.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to prevent contamination.
- Use the QuikClot products as directed to control bleeding effectively.
- Ensure that the Emergency Trauma Dressing is applied correctly to avoid further injury.
- The hypothermia blanket should be used to maintain body temperature in emergency situations.
- Compressed gauze should be used to pack wounds as needed.
- The Nasopharyngeal Airway and Decompression Needle are for trained personnel only; improper use may cause harm.

Instructions for Installation and Usage

1. Accessing the Kit:

- Open the MOLLE pouch by pulling the top flap.
- Ensure that the contents are organized and easily accessible.

2. Using Trauma Gloves:

- Roll out the gloves from the pouch.
- Ensure your hands are dry before putting on the gloves.
- Dispose of gloves properly after use.

3. Applying QuikClot Bleeding Control Dressing:

- Remove the dressing from its packaging.
- Apply directly to the bleeding area and hold in place for at least 3 minutes.
- Monitor for continued bleeding; if necessary, apply additional dressings.

4. Using the Emergency Trauma Dressing:

- Open the dressing and place it over the wound.
- Secure it in place using the attached ties.
- Do not remove if the dressing becomes soaked; apply additional dressings on top.

5. Using the Hypothermia Blanket:

- Unfold the blanket and wrap it around the person in need.
- Ensure the shiny side is facing inward to retain body heat.

6. Using Compressed Gauze:

- Remove the gauze from its packaging.
- Pack it into the wound as necessary to control bleeding.

7. Using the Nasopharyngeal Airway:

- Only use if trained to do so.
- Lubricate the airway before insertion and gently insert it into the nostril.

8. Using the Decompression Needle:

- Only use if trained to do so.
- Ensure proper technique to avoid injury; consult a medical professional if unsure.

Disposal Instructions

- Dispose of used gloves, dressings, and gauze in accordance with local regulations for biohazard waste.
- Check local guidelines for proper disposal of medical supplies and products.
- Recycle packaging materials where possible.

Contact Information for Further Support

For any questions or concerns regarding the TRAUMA KIT NOW! SMALL, please consult the manufacturer's website or customer service for assistance. Ensure that you keep this guide for future reference.

Thank you for prioritizing safety with the TRAUMA KIT NOW! SMALL.

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