TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES-MULTICAM

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES-MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101023
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.44kg
- Shipping height: 51mm
- Shipping width: 203mm
- Shipping length: 279mm
- UPC: 810073653283

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Trauma Kit NOW! Small Trauma Kit Safety Instructions

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL by Blue Force Gear. This kit is designed to provide essential medical supplies in a compact and efficient manner. Please read this safety instruction guide carefully to ensure safe and effective use of the kit.

General Safety Guidelines

- **Product Safety:** This kit is designed for safe use. Ensure that you familiarize yourself with the contents and their purposes.
- **Enhanced Recalls:** In the unlikely event of a recall, you will be notified through standardized notices. Free remedies will be provided for any dangerous products.
- **Online Shopping:** If purchased online, you have the same safety protections as if you bought it in a physical store.
- **Special Consumer Focus:** This kit is not intended for children. Keep out of reach of minors.
- **EU Contact Point:** For safety inquiries, please refer to the appropriate EU-based contact.
- **Rapid Alerts:** Stay informed about any unsafe products through the EU's Safety Gate system.

Specific Safety Precautions for Use

- **Avoid Contact with Eyes:** Ensure that all products, especially those containing chemicals like QuikClot, are kept away from the eyes.
- **Proper Use of Gloves:** Always use the trauma gloves when handling any medical supplies to maintain hygiene and prevent contamination.
- **Check Expiration Dates:** Regularly check the expiration dates on all medical supplies and replace any expired items.
- **Follow Instructions:** Use each item as directed in its respective instructions to avoid misuse and potential harm.
- **Store in a Safe Place:** Keep the kit in a location that is easily accessible but safe from children and pets.

Instructions for Installation and Usage

- **Initial Setup:**
 - Open the kit and familiarize yourself with all included items.
 - Ensure that the kit is stored in a dry and cool place to maintain the integrity of the supplies.
- **Using the Supplies:**
 - **Trauma Gloves:** Roll out the gloves and wear them before tending to any injuries.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound and apply pressure as needed.
 - **Emergency Trauma Dressing:** Use to cover larger wounds; follow the instructions on the packaging for application.
 - **Hypothermia Blanket:** Unfold and wrap around the injured individual to maintain body heat.
 - **Flat Fold Tape:** Use to secure dressings in place.
 - **Compressed Gauze:** Apply to wounds and secure with tape or other means.
 - **Chest Seal:** Follow the packaging instructions for application on chest wounds.
 - **Nasopharyngeal Airway:** Insert carefully as per instructions, ensuring it is used only by trained individuals.
 - **Decompression Needle:** Use only in emergencies and by trained personnel.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations.
- Do not throw away medical supplies in regular trash; check for local hazardous waste disposal options.
- Ensure that all items are disposed of safely to prevent accidental use or exposure.

Contact Information for Further Support

For any questions or concerns regarding the Trauma Kit NOW! SMALL, please refer to the provided contact information on the packaging or visit the manufacturer's website for assistance.

Thank you for prioritizing safety. Proper use of this kit can make a significant difference in emergency situations.

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