TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL-MOLLE-PRO SUPPLIES-RANGER GREEN

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL-MOLLE-PRO SUPPLIES-RANGER GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101024
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.544kg
- Shipping height: 51mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653290

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Trauma Kit NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL from Blue Force Gear. This compact medical kit is designed for quick access to essential trauma supplies. It is important to follow safety guidelines to ensure effective use and to minimize risks during emergencies.

General Safety Guidelines

- Always check the contents of your trauma kit regularly to ensure all items are present and in good condition.
- Store the kit in a cool, dry place away from direct sunlight.
- Keep the kit out of reach of children unless under adult supervision.
- Familiarize yourself with the contents and their uses before an emergency occurs.
- Ensure that all users are trained in basic first aid and CPR.
- Report any unsafe products or accidents to the appropriate authorities.
- Check for recall updates on the EU's Safety Gate platform to stay informed about potential hazards.

Specific Safety Precautions for Use

- **Trauma Gloves:** Use gloves to prevent contamination and protect yourself and the patient from infection.
- **QuikClot Products:** Follow the instructions provided with the QuikClot products carefully to ensure proper application for bleeding control.
- **Emergency Trauma Dressing:** Ensure the dressing is applied correctly to avoid further injury; do not apply directly to open wounds without proper training.
- **Hypothermia Blanket:** Use the blanket to maintain body temperature, especially in cases of shock or exposure.
- **Compressed Gauze:** Use gauze to control bleeding; apply firmly and securely.
- **Chest Seal:** Ensure the chest seal is applied correctly to prevent air from entering the chest cavity.
- **Nasopharyngeal Airway and Decompression Needle:** Only use these items if you are trained to do so, as improper use can cause harm.

Instructions for Installation and Usage

1. **Open the Kit:** Unzip the pouch and lay out all items for easy access.
2. **Assess the Situation:** Determine the nature of the injury and the necessary supplies.
3. **Put on Trauma Gloves:** Always wear gloves before providing care.
4. **Control Bleeding:** Use the QuikClot product or compressed gauze to control any bleeding.
5. **Apply Dressings:** Use the emergency trauma dressing as needed.
6. **Maintain Body Temperature:** Wrap the patient in the hypothermia blanket if necessary.
7. **Monitor the Patient:** Keep an eye on the patient's condition and be ready to provide further assistance.
8. **Use the Nasopharyngeal Airway:** Only if trained, insert the airway to maintain an open airway.
9. **Decompression Needle:** Only use if trained and in life-threatening situations to relieve tension pneumothorax.
10. **Clean Up:** Dispose of used materials properly and sanitize your hands after use.

Disposal Instructions

- Dispose of used gloves, dressings, and other materials in accordance with local regulations for medical waste.
- For nonhazardous items, follow standard waste disposal guidelines.
- Always ensure that any sharps (e.g., needles) are disposed of in a designated sharps container.

Contact Information for Further Support

For questions or concerns regarding the Trauma Kit NOW! SMALL, please reach out to your local distributor or refer to the product packaging for further information.

This safety instruction guide is designed to help you use the Trauma Kit NOW! SMALL safely and effectively. Always prioritize safety and training to ensure the best outcomes in emergency situations.

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