TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL-MOLLE- PRO SUPPLIES-WOLF GRAY

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL-MOLLE- PRO SUPPLIES-WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101025
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.442kg
- Shipping height: 102mm
- Shipping width: 178mm
- Shipping length: 254mm
- UPC: 810073653306

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TRAUMA KIT NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL, the latest medical kit from Blue Force Gear. This compact and versatile kit is designed for effective emergency response, offering essential supplies in a MOLLEcompatible pouch. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace items as necessary.
- Keep the kit out of reach of children and vulnerable individuals unless supervised by an adult.
- Familiarize yourself with the contents and their uses before an emergency arises.
- In case of an emergency, prioritize calling for professional medical assistance.

Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear gloves when providing first aid to prevent contamination.
- **QuikClot Dressing:** Apply pressure to control bleeding, but do not use if the wound is deep or if there is a foreign object embedded.
- **Emergency Trauma Dressing:** Ensure the dressing is applied securely to promote effective wound closure.
- **Hypothermia Blanket:** Use the blanket to prevent heat loss in injured individuals, especially in cold environments.
- **Nasopharyngeal Airway:** Only use if trained to do so; improper insertion can cause injury.
- **Decompression Needle:** Use only if trained; improper use can lead to severe injury or complications.

Instructions for Installation and Usage

1. **Familiarization:** Before use, review the contents of the kit and understand the purpose of each item.
2. **Accessing the Kit:** Open the MOLLE pouch by pulling the top flap. The singlemotion design allows for quick access.
3. **Using Supplies:**
 - **Trauma Gloves:** Roll out two pairs and wear them before handling any medical supplies.
 - **Bleeding Control:** Use the QuikClot dressing by placing it directly on the wound and applying firm pressure.
 - **Dressing Application:** For the Emergency Trauma Dressing, unwrap and place it over the wound, securing it with the provided ties.
 - **Heat Retention:** Wrap the Hypothermia Blanket around the individual to maintain body temperature.
 - **Airway Management:** If needed, carefully insert the Nasopharyngeal Airway as per training guidelines.
 - **Needle Decompression:** If trained, use the Decompression Needle in cases of suspected tension pneumothorax.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations for medical waste.
- Do not attempt to reuse singleuse items such as gloves, dressings, or gauze.
- Ensure that all waste is sealed in a suitable container to prevent contamination.

Contact Information for Further Support

For any inquiries regarding the TRAUMA KIT NOW! SMALL, please refer to the EU contact point listed on the packaging or product documentation.

Conclusion

The TRAUMA KIT NOW! SMALL is designed to provide essential medical supplies in a compact form. By following these safety instructions, you can ensure that you and others can respond effectively in emergencies. Always prioritize safety and seek professional medical help when necessary.

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