TRAUMA KIT NOW! - SMALL - TRAUMA KIT NOW! SMALL - MOLLE - ADVANCED SUPPLIES-MULTICAM

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: TRAUMA KIT NOW! SMALL - MOLLE - ADVANCED SUPPLIES-MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101028
- Mfr. No.: HW-TKN-SMALL-AD
- Size: 1 Person
- Delivery weight: 0.567kg
- Shipping height: 279mm
- Shipping width: 76mm
- Shipping length: 229mm
- UPC: 810073653337

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TRAUMA KIT NOW! SMALL Safety Instructions Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL by Blue Force Gear. This advanced medical pouch is designed to provide essential supplies for emergency situations. It is crucial to familiarize yourself with the safety guidelines and instructions to ensure safe and effective use of the kit.

General Safety Guidelines

- Ensure the kit is stored in a location that is easily accessible in case of an emergency.
- Regularly check the contents of the kit to ensure all items are present and within their expiration dates.
- Keep the kit out of reach of children to prevent accidental use.
- If you are unsure about how to use any item in the kit, seek guidance from a qualified medical professional.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear Trauma Gloves when handling blood or bodily fluids to minimize the risk of infection.
- Use the QuikClot products as directed in the instructions provided with the product to ensure proper use.
- Be cautious when using sharp objects, such as scissors or needles, included in the kit.
- Ensure that the Emergency Trauma Dressing is applied correctly to avoid further injury.
- If using the Nasopharyngeal Airway, ensure it is the correct size for the individual and is inserted gently to avoid injury.

Instructions for Installation and Usage

1. Opening the Kit

- To access the contents, pull out the top of the pouch using a single motion.
- Ensure you have a clear area to work in when opening the kit.

2. Using Trauma Gloves

- Roll out two pairs of Trauma Gloves from their packaging.
- Ensure your hands are clean before putting on the gloves.

3. Applying QuikClot Bleeding Control Dressing

- Open the QuikClot package and apply it directly to the bleeding area.
- Apply pressure to the dressing for at least 3 minutes or until bleeding stops.

4. Using the Emergency Trauma Dressing

- Unwrap the dressing and apply it to the wound.
- Secure it in place using the attached ties.

5. Using the Hypothermia Blanket

- Remove the blanket from its packaging and unfold it.
- Wrap the patient securely to maintain body temperature.

6. Using Flat Fold Tape

- Cut the desired length of tape and apply it to secure dressings or gauze.
- Ensure the tape does not restrict blood flow.

7. Using Compressed Gauze

- Unwrap the gauze and pack it into the wound if necessary.
- Apply pressure to control bleeding.

8. Using the Nasopharyngeal Airway

- Select the appropriate size for the individual.
- Lubricate the airway and gently insert it into the nostril, following the natural curve.

9. Using the Decompression Needle

- Ensure the area is clean and identify the site for insertion.
- Follow specific instructions provided with the needle for safe use.

Disposal Instructions

- Dispose of used gloves, dressings, and other single-use items in a sealed bag to prevent contamination.
- Follow local regulations for the disposal of medical waste.
- Ensure that any sharp objects are disposed of in a designated sharps container.

Contact Information for Further Support

For any questions regarding the TRAUMA KIT NOW! SMALL, please refer to the product manual or contact the manufacturer's customer support. Always ensure that you are following the latest safety guidelines and product instructions.

By adhering to these safety guidelines and instructions, you can effectively use the TRAUMA KIT NOW! SMALL to assist in emergency situations while ensuring your safety and the safety of others.

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